

Assay	\$/sample	Company	Sample size	Sample	Sensitivity	Assay Flavor	%CV	
							intra	inter
A1C	\$7.88	Med. Computer Sys.	50 ul	Blood	<1%CV	Analyzer		2.34
ACTH	\$11.71	DSL	200 ul in dup.	EDTA plasma	1.3pg/ml	IRMA		
Adiponectin ELISA	\$12.50	Linco	100 ul	Serum/Plasma	0.78ng/ml	Elisa		
Adiponectin HMW Elisa	\$19.44	Linco	50 ul	Serum/Plasma		Elisa		
Adiponectin RIA	\$8.50	Linco	100 ul	Serum/Plasma	0.5ng/ml	RIA		
Adrenaline	\$12.66	ALPCO	1.0ml	EDTA plasma	0.9ng/ml	RIA	-	3.17
Androstenedione EIA	\$7.50	DSL	50 ul in dup.	Serum/Plasma	0.03ng/ml	RIA		
CICP(PICP)	\$18.57	QUIDEL	50 ul in dup.	Serum		Elisa	6.08	
CMVlgG Elisa	\$7.64	United Biotech	50 ul	Plasma/serum		Elisa		
Cortisol serum	\$4.29	DSL	25 ul in dup.	Serum/Plasma	0.11ug/dl	RIA	4.65	6.3
Cortisol salivary	\$4.00	DSL	25 ul in dup.	Salivary	0.11ug/dl	Elisa	8.22	5.03
Cortisol urine	\$4.29	DSL	250 ul	Urine	0.11ug/dl	RIA	4.61	8.01
C-Peptide	\$2.30	Linco	100 ul in dup.	Plasma/serum	0.1ng/ml	RIA	2.47	0.61
Creatinine serum	\$2.25	Beckman	25 ul in dup.	Serum	L.R.0.2-25mg/dl	Analyzer	0.47	3.1
Creatinine urine	\$2.25	Beckman	25 ul in dup.	Urine	L.R.0.2-25mg/dl	Analyzer	0	1.7
CrossLaps(serum)	\$22.86	Nordic Bioscience Diagnostics A/S	50 ul in dup.	Serum	95% cof.interval	Elisa	8.15	0.01
CRP	\$10.95	ALPCO	50 ul	Plasma/serum	0.5 ng/ml	Elisa	6.33	2.2
DHEA	\$7.50	DSL	100 ul in dup.	Plasma/serum	0.02ng/ml	RIA	4.31	5.53
DHEA-S	\$4.89	ALPCO	25 ul in dup.	Plasma/serum	0.02ug/ml	Elisa	3.28	5.78
D-Pyridinoline	\$20.83	QUIDEL	50 ul in dup.	Urine	1 nM	Elisa	8.44	7.79

Assay	\$/sample	Company	Sample size	Sample	Sensitivity	Assay Flavor	%CV	
							intra	inter
Erythropoietin	\$19.86	R&D Systems	100 ul in dup.	Plasma/serum	0.6mIU/ml	Elisa	7.56	4.53
Estradiol(UL.sen.)	\$4.72	DSL	200 ul in dup.	Plasma/serum	2.2pg/ml	RIA	4.37	7.54
Estrone DA	\$8.61	DSL	50 ul in dup.	Plasma/serum	1.2pg/ml	RIA	5.45	6.71
FGF-23(Intact)	\$17.64	Alpco	150 ul in dup.	EDTA plasma	recovery 90-110%	Elisa		
FSH	\$3.89	DSL	100 ul in dup.	Serum	0.2 mIU/ml	IRMA	3.57	5.76
GHRELIN Total	\$9.19	Linco	100 ul in dup.	Serum/Plasma	100 pg/ml	RIA	5.17	4.76
Glucose	\$1.25	Beckman	10 ul in dup.	Serum	L.R.1-450mg/dl	Analyzer	1.38	1.52
Helicopeptide(urine)	\$21.43	QUIDEL		Urine		Elisa	3.69	6.43
HGH	\$7.50	Nichols	100 ul in dup.	Serum	0.02 ng/ml	RIA	8.31	0.76
Homocysteine	\$15.29	ALPCO	50 ul	EDTA plasma	1.5umol/L	Elisa		
ICTP	\$15.28	DiaSorin	100 ul in dup.	Serum	0.5ug/L	RIA	7.84	2.5
IFN -Y	\$12.29	R&D Systems				Elisa		
IGF-I	\$14.00	Alpco	50 ul	Plasma/serum	0.1ng/ml	RIA	7.15	9.64
IL-1 ra	\$14.14	R&D Systems	200 ul in dup.	Serum/Plasma	22pg/ml	Elisa		
IL-10	\$16.14	R&D Systems	200 ul in dup.	Serum/Plasma	<2 pg/ml	Elisa		
IL-10 HS	\$17.00	R&D Systems	200 ul in dup.	Serum/Plasma	<0.5 pg/ml	Elisa		
IL-11	\$13.29	R&D Systems	100 ul in dup.	Serum/Plasma	8.0pg/ml	Elisa		
IL-12-bioactive	\$13.86	R&D Systems	200 ul in dup.	Serum/Plasma	<5 pg/ml	Elisa		
IL-12-inhibitor antagonest	\$13.86	R&D Systems	100 ul in dup.	Serum/Plasma	<15 pg/ml	Elisa		
IL-13	\$12.29	R&D Systems	100 ul in dup.	Serum/Plasma	<32 pg/ml	Elisa		

Assay	\$/sample	Company	Sample size	Sample	Sensitivity	Assay Flavor	%CV	
							intra	inter
IL-1HS Quantikine	\$17.00	R&D Systems	150 ul in dup.	Serum/Plasma	0.1pg/ml	Elisa		8.5
IL-4 HS	\$16.00	R&D Systems	200 ul in dup.	Serum		Elisa		
IL-6 HS Quantikine	\$17.00	R&D Systems	100 ul in dup.	Serum/Plasma	0.039pg/ml	Elisa	5.13	4.97
IL-6sR	\$14.14	R&D Systems	50 ul	Serum/Plasma	1.5-15.1pg/ml	Elisa	4.43	3.54
IL-8	\$14.00	R&D Systems	50 ul in dup.	Serum/Plasma	3.5pg/ml	Elisa		
INF-gamma	\$9.09	R&D	50 ul in dup.	Plasma/serum	8 pg/mL	Elisa		
Insulin	\$2.19	Linco	100 ul in dup.	Serum/Plasma	2 uU/ml	RIA	4.42	5.95
IPTH	\$10.28	ALPCO	25 ul in dup.	Serum/EDTA Plasma	0.9pg/ml	Elisa	4.09	4.99
Leptin	\$3.56	Linco	100 ul in dup.	Plasma/serum	0.5 ng/ml	RIA	3.85	3.81
LH	\$3.89	DSL	100 ul in dup.	Serum	0.1mIU/ml	IRMA	3.37	5.9
LH Active Elisa	\$6.94	DSL				Elisa		
Lipid profile	\$10.25	Med. Computer Sys.	35 ul	Serum	0.8%	Analyzer	3.082	3.72
LP(a)	\$25.57	Diapharma	50 ul	serum/citrate plasma	0.001mg/dl	Elisa	3.69	9.58
MCP-1	\$13.14	R&D Systems	250 ul	Plasma/serum	5.0pg/ml	Elisa	8.67	7.01
Metra BAP	\$20.83	Quidel	20 ul in dup.	Serum	0.7U/L	Elisa	5.49	5.83
MIF	\$13.75	R&D Systems				Elisa		
MMP-9(total)	\$13.14	R&D Systems	10 ul	Serum/Heparin Plasma	0.156ng/ml	Elisa	3.91	8.99
Myeloperoxidase(MPO)	\$16.39	R&D Systems	100ul pasma(dil 1:30)	Plasma	1.5ng/ml	Elisa		
NEFAC	\$3.67	Wako Chem.	25 ul	Serum	L.R.0.25-1mEo/L	in house	5.73	4.99
Nitrotyrosine HBT	\$7.41	Cell Sciences	50 ul	Serum/Plasma		Elisa		

Assay	\$/sample	Company	Sample size	Sample	Sensitivity	Assay Flavor	%CV	
							intra	inter
Noradrenaline	\$12.66	ALPCO	1.0ml	EDTA plasma	4.5ng/ml	RIA	-	4.86
Noradrenaline-serum	\$18.29	ALPCO	300ul	EDTA plasma	4.5ng/ml	RIA		7.73
NTX	\$15.00	OSTEX	25 ul in dup.	Urine	Recovery98%	Elisa	4.56	8.83
NTXS	\$35.17	OSTEX	50 ul	Serum	Recovery98%	Elisa	1.1	0.83
OPG	\$21.53	ALPCO	50 ul in dup.	Plasma/serum		Elisa		
Osteocalcin	\$10.50	Immutopics	10 ul in dup.	Serum	0.06 ng/ml	IRMA	2.97	6.3
PAI/tPA Complex	\$12.80	Alexis Corp.	50 ul	Plasma/serum	1-15ng/ml	Elisa	5.67	0.84
PAI-1	\$15.00	Alexis Corp.	25 ul in dup.	EDTA or citrate plasma	0.5 ng/ml	Elisa		
PIIINP	\$16.67	DiaSorin	200 ul in dup.	Serum is recom./EDTA Plasma		RIA		
PINP	\$16.67	IDS Inc.	50 ul in dup.	Serum	2.0ug/L	RIA	4.14	2.74
Progesterone	\$4.72	DSL	25 ul in dup.	Plasma/serum	0.05ng/ml	RIA	7.3	8.96
Proinsulin	\$2.76	Linco	200 ul in dup.	Plasma/serum	2.0pM	RIA	2.65	6.62
Prolactin	\$3.89	DSL	25 ul in dup.	Serum	0.1ng/ml	RIA	2.24	1.52
Prolactin Active Elisa	\$6.94							
Pro-mmp-1	\$12.86	R&D Systems	100 ul in dup.	Plasma/serum	0.021ng/ml	Elisa		
PYD(serum)	\$20.83	QUIDEL		Serum		Elisa		
PYD(urine)	\$20.83	QUIDEL		Urine		Elisa		
Resistin	\$14.63	Linco	20 ul in dup.	Serum/EDTA Plasma	0.16ng/ml	Elisa	3.67	3.16
serum Amyloid A	\$11.94	Anigenix America		Serum	<1.2ng/ml	Elisa		
sE-Selectin	\$13.57	R&D Systems	25 ul	Serum/Heparin Plasma	0.1 ng/ml	Elisa		

Assay	\$/sample	Company	Sample size	Sample	Sensitivity	Assay Flavor	%CV	
							intra	inter
sgp130	\$14.14	R&D Systems	50 ul	Serum/Plasma	0.08ng/ml	Elisa	6.8	1.08
SHBG Active Elisa	\$12.50	DSL		Serum/Plasma		Elisa		
SHBG Active IRMA	\$12.50	DSL	25 ul in dup.	Serum/Plasma	3.0nmol/L	IRMA	6.54	7.56
sICAM-1	\$13.29	R&D Systems	15 ul	Serum/Plasma	0.35ng/ml	Elisa	2.48	1.83
sP-Selectin	\$13.29	R&D Systems	25 ul	Serum/Heparin Plasma	0.5ng/ml	Elisa	3.64	6.87
sRANKL	\$27.43	ALPCO	100 ul in dup.	Plasma/serum		Elisa		
sTNF receptor 2	\$14.14	R&D Systems	100 ul	Serum/Plasma	0.6pg/ml	Elisa	2.55	5.93
sTNFreceptor 1	\$14.14	R&D Systems	100 ul	Serum/Plasma	0.77pg/ml	Elisa	1.31	8.95
sV-CAM-1	\$13.75	R&D Systems	50 ul	Serum/Plasma	2.0ng/ml	Elisa	1.89	5.97
T4	\$3.06	DSL	50 ul in dup.	Serum/Plasma	0.08ng/dL	RIA		
Testosterone Free	\$6.39	DSL	50 ul in dup.	Serum	0.18pg/ml	RIA		
Testosterone Total	\$4.72	DSL	50 ul in dup.	Serum/Plasma	0.08ng/ml	RIA		
TGF beta 1	\$14.14	R&D Systems	200 ul in dup.	Serum/Plasma		Elisa		
TIMP-1	\$12.50	R&D Systems	25ul	Serum/Plasma	0.08ng/ml	Elisa		
TNF-alpha	\$17.00	R&D	200 ul in dup.	Plasma/serum	0.06-0.32pg/ml	Elisa	8.22	6.05
tPA	\$13.11	Alexis Corp.	25 ul in dup.	EDTA or citrate plasma		Elisa		9.92
TSH	\$6.11	DSL	100 ul in dup.	Serum	0.4mIU/ml	IRMA	2.1	5.2
Ultra sen. GH	\$5.56	DSL	100 ul in dup.	Serum	0.7pg/ml	RIA	8.41	5.01
Ultra sen. insulin	\$2.19	Linco	50 ul in dup.	Serum/Plasma	0.2 uU/ml	RIA	2.69	3.27
Vit D	\$13.19	DiaSorin	50 ul	Serum/Plasma	3.0 ng/ml	RIA	5.19	7.9

Assay	\$/sample	Company	Sample size	Sample	Sensitivity	Assay Flavor	%CV	
							intra	inter
vWF Activity	\$12.86	American diagnostica inc.	100 ul in dup.	EDTA or citrate plasma		Elisa	2.57	8.82
YKL 40	\$21.43	QUIDEL	20 ul in dup.	Serum	20ng/ml	Elisa	4.22	-